

# gloveon Matador



## Nitrile Exam Gloves Powder Free, Long Cuff

With quality, consistency and reliability built into the GloveOn's 8 mil nitrile glove, even a bull fighter will have the courage, trust and confidence to go into the ring to face the challenges while donning the **Matador**.



| GloveOn Matador                    |                         |                     |
|------------------------------------|-------------------------|---------------------|
| <b>Length (mm)</b>                 | ≥ 280                   |                     |
| <b>Thickness Measurements (mm)</b> |                         |                     |
| Palm (centre of Palm)              | 0.20 ± 0.02             |                     |
| Finger (13mm ± 3mm from tip)       | 0.24 ± 0.02             |                     |
| <b>Physical Properties</b>         | <b>Before Ageing</b>    | <b>After Ageing</b> |
| Tensile Strength (MPa)             | ≥ 18                    | ≥ 16                |
| Elongation (%)                     | ≥ 500                   | ≥ 400               |
| <b>Inspection Levels &amp; AQL</b> | <b>Inspection Level</b> | <b>AQL</b>          |
| Watertightness                     | G1                      | 1.50                |
| Physical Dimensions                | S2                      | 4.00                |
| Tensile Strength                   | S2                      | 4.00                |
| Visual Inspection (Major)          | S4                      | 2.50                |
| Visual Inspection (Minor)          | S4                      | 4.00                |
| Particulate Residue                | N = 5                   | ≤ 2mg/glove         |

Tested in accordance with EN 374-3:2003 with the following chemical permeation performance levels:

| Chemical                    | Performance Level | Permeation Performance Level | Measured Breakthrough Time (Minutes) |
|-----------------------------|-------------------|------------------------------|--------------------------------------|
| 4% Chlorhexidine gluconate  | 6                 |                              |                                      |
| 40% Sodium hydroxide        | 6                 |                              |                                      |
| 35% Gluteraldehyde          | 6                 | 1                            | > 10                                 |
| 1% Gluteraldehyde           | 6                 | 2                            | > 30                                 |
| 4 % Gluteraldehyde          | 6                 | 3                            | > 60                                 |
| 10% Formalin (Formaldehyde) | 6                 | 4                            | > 120                                |
| 35% Ethanol                 | 3                 | 5                            | > 240                                |
| n-Hexane                    | 2                 | 6                            | > 480                                |
| Isopropanol                 | 3                 |                              |                                      |
| Acetone                     | Not recommended   |                              |                                      |
| 3% Hydrogen Peroxide        | 6                 |                              |                                      |
| n-Heptane (J)               | 5                 |                              |                                      |
| 96% Sulphuric Acid (K)      | 2                 |                              |                                      |
| 40% Sodium Hydroxide (L)    | 6                 |                              |                                      |

### REORDER CODE

|         |         |
|---------|---------|
| NTR18XS | X-SMALL |
| NTR18SS | SMALL   |
| NTR18MM | MEDIUM  |
| NTR18LL | LARGE   |
| NTR18XL | X-LARGE |

### FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Long cuff
- Sky blue colour

### PACKAGING

50 gloves per box for XS to XL  
10 boxes per carton

### REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k),  
MDD 93/42/EEC, REACH,  
ROHS Directive 2002/95/EC,  
PPE 89/686/EEC

### STANDARDS

ASTM D6319, EN 455 part 1, 2, 3 & 4,  
EN 374 part 1, 2, 3, AS/NZS 4011:1997

### MANUFACTURING ACCREDITATIONS

ISO 9001:2008  
ISO 13485:2003  
EN ISO 13485:2003

A brand by

**D:mun**  
www.munglobal.com

**Mun Australia**  
Toll free: 1800 456 837 (GLOVES)  
info.au@munglobal.com

  
Proud Supporter of  
Birthing Kit Foundation (Australia)

 Protection Always On